

Coefficient of Friction Test Fixture


TH62-M and TH62-F



Fixture to test friction coefficient
Pulley with light moving ball bearing


Tensile force: 500 N
Platen: 270x150 mm
470x150 mm
570x150 mm

TH62-M:
ASTM-D1894, TAPPI T549




Fixture including weight
(covered with cellular rubber)
200 g

TH62-F:
DIN53375, ISO8295




Fixture including weight
(covered with felt)
200 g



Alu 3 kg

TH62-M




TH62-F

Optional Accessory for TH62



TH62-HS

Folding sledge with magnet fastener to keep thin samples (like foils, paper) secure in place



TH62-HS

Similar Fixtures Conform to Other Standards




THS62-20-500


Peel cling of stretch wrap film fixture
ASTM-D5458
Size 500x150 mm

THS62-20-600

Peel cling of stretch wrap film fixture
ASTM-D5458
Size 600x150 mm



TH62-20-500



TH62-20-600



TH62-40-ISO21182

Slide weight, 50 N