

Coefficient of Friction Test Fixture


TH62-M and TH62-F



Fixture to test friction coefficient
Pulley with light moving ball bearing


Tensile force: 500 N
Platen: 270x150 mm
470x150 mm
570x150 mm

TH62-M:
ASTM-D1894, TAPPI T549




Fixture including weight
(covered with cellular rubber)
200 g

TH62-F:
DIN53375, ISO8295



Fixture including weight
(covered with felt)
200 g




Alu 3 kg

TH62-M

AMSE S.r.l. Analytical & Material Science Equipment

P. IVA 09762820018 | Strada della Cebrosa 86 | 10156 Torino (TO) Italy | Tel. +39.011.27.30.000 | Fax +39.011.27.404.58 |

info@amse.it www.amse.it




TH62-F

Optional Accessory for TH62



TH62-HS

Folding sledge with magnet fastener to keep thin samples (like foils, paper) secure in place



TH62-HS

AMSE S.r.l. Analytical & Material Science Equipment

P. IVA 09762820018 | Strada della Cebrosa 86 | 10156 Torino (TO) Italy | Tel. +39.011.27.30.000 | Fax +39.011.27.404.58 |

info@amse.it www.amse.it

Similar Fixtures Conform to Other Standards




THS62-20-500


Peel cling of stretch wrap film fixture
ASTM-D5458
Size 500x150 mm

THS62-20-600

Peel cling of stretch wrap film fixture
ASTM-D5458
Size 600x150 mm



TH62-20-500



TH62-20-600



TH62-40-ISO21182

Slide weight, 50 N

AMSE S.r.l. Analytical & Material Science Equipment

P. IVA 09762820018 | Strada della Cebrosa 86 | 10156 Torino (TO) Italy | Tel. +39.011.27.30.000 | Fax +39.011.27.404.58 |

info@amse.it www.amse.it